

LED DAYLIGHT LAMP



**Slim compact
design**



**USB power
supply**

**Adjustable
brightness
and light
temperature**



The LED daylight lamp imitates natural sunlight and provides a dose of daylight in the dark winter months. It is the perfect companion for the home or in the office.

www.drivedevilbiss.co.uk

drive

DeVilbiss
HEALTHCARE

LED daylight lamps produce a light intensity similar to sunlight and can be used to supplement missing exposure during the winter months.

Light boxes come in a variety of designs that produce a very bright light. The intensity of the light is measured in lux.

Most people can use light therapy safely. Light boxes have filters that remove harmful ultraviolet (UV) rays, so there's no risk of skin or eye damage for most people.



**Up to
10,000
Lux light
intensity**

Features:

- » LED Technology
- » Adjustable brightness from 2,000 lux to 10,000 lux
- » Adjustable light temperature from 3000K to 6500K (bright white to warm yellow)
- » Touch control for hue dimming
- » Up to 10,000 Lux light intensity (10cm)
- » Compact and portable
- » Fold away stand
- » Automatic timing setting of 15, 30, 45, 60 minutes
- » USB power supply



**Compact
and portable
with foldaway
stand**

Specifications:



Product Code:
PCWB001

Height:
24cm

Width:
14.6cm

Depth:
1.7cm

Product Weight:
658g



www.drivedevilbiss.co.uk

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.